

CT Distancing



Delivering
Mobile CT
Solutions
Since 2005

Keep Your Most Vulnerable Patients Safe with OmniTom® Mobile Imaging at the Bedside

In today's challenging healthcare environment, OmniTom brings the power of state-of-the-art, multi-slice CT right to the bedside to safeguard your most vulnerable patients.

Minimize the exposure of pediatric patients to staff and visitors by eliminating transport to the CT Suite. Keep PICU and other fragile patients attached to their vital life-sustaining equipment throughout imaging procedures.

- High quality 16-slice images
- 40 cm bore for versatile imaging of small patients
- Self-shielded design
- Battery-operation
- Immediate image viewing
- Easily maneuverable from room to room



With its dedicated scanning platforms, OmniTom can image neonatal and pediatric patients.

The mobile scanner is an indispensable tool for collecting real time data on critically ill patients as well as for infectious disease screening, back up scanning and a full range of point-of-care needs.

Bring the power of Mobile CT to your patients.

For more information, email: MobileCT@NeuroLogica.com or call: 978-716-2726

Our mobile CT systems are proudly manufactured in the USA

© 2020 NeuroLogica Corporation. All rights reserved. Danvers, MA